

# SIP goes Proteomics – Elucidation of structure and function of microbial communities

PD Dr. Martin von Bergen

Keystone 2010



Microbial Life is complex

→ Who eats what, where and when?

Courtesy of Matthias Kästner, UFZ

12C-Toluol + Nitrat Aussenseite 11.06.07 1kV 5mm

2 μm

# History and features of SIP approaches in microbiology

1998: Labelling of phospholipids (Boschker et al.)

- Direct linking of microbial populations to specific biogeochemical processes by  $^{13}\text{C}$ -labelling of biomarkers
- Acetate consumption rather by SRB similar to *Desulfotomaculum acetoxidans* than by *Desulfobacter* spp.
- + Combined assessment of species and metabolic function
- Limited taxonomic resolution

# History and features of SIP approaches in microbiology

1998: Labelling of phospholipids (Boschker et al.)

- + Combined assessment of species and metabolic function
- Limited taxonomic resolution

2000: DNA-SIP (Radajewski et al.)

- Labelled DNA can be purified, amplified and sequenced yielding wide taxonomic information coupled with information on metabolic activity
- Example of different methylotrophic species obtained by feeding  $^{13}\text{CHOH}$
- + Combination of genetic methods allows even culture independent assessment of species
- Information on the potential but not the actual function
- Density gradient centrifugation enables only low resolution, results in high minimal labelling

# History and features of SIP approaches in microbiology

1998: Labelling of phospholipids (Boschker et al.)

- + Combined assessment of species and metabolic function
- Limited taxonomic resolution

2000: DNA-SIP (Radajewski et al.)

- + Combination of genetic methods allows even culture independent assessment of species
- Information on the potential but not the actual function
- Density gradient centrifugation provides only low resolution, results in high minimal labelling

2002: RNA-SIP (Manefield et al.)

- Labelled RNA can be purified, amplified and sequenced yielding wide taxonomic information coupled with information on functional aspects
- Example of phenol degradation in aerobic bioreactor and identification of a *Thauera* spp. as a keyplayer
- + Combination of genetic methods allows even culture independent assessment of species
- + Detection of transcription of functional genes
- Density gradient centrifugation enables only low resolution, results in high minimal labelling

# History and features of SIP approaches in microbiology

1998: Labelling of phospholipids (Boschker et al.)

- + Combined assessment of species and metabolic function
- Limited taxonomic resolution

2000: DNA-SIP (Radajewski et al.)

- + Combination of genetic methods allows even culture independent assessment of species
- Information on the potential but not the actual function
- Density gradient centrifugation provides only low resolution, results in high minimal labelling

2002: RNA-SIP (Manefield et al.)

- + Combination of genetic methods allows even culture independent assessment of species
- + Detection of transcription of functional genes
- Density gradient centrifugation enables only low resolution, results in high minimal labelling

2006: Nano-SIMS (Lechene et al.)

- High spatial resolution in combination with measurement of secondary ions
- Analysis of microbial nitrogen fixation by detection of  $^{15}\text{N}$  incorporation
- + Single cell analysis
- + high dynamic range
- + Direct correlation between metabolomic activity and species identification in combination with FISH
- Due to tedious sample preparation only consortia of limited complexity can be analysed
- relatively costly instruments

# History and features of SIP approaches in microbiology

1998: Labelling of phospholipids (Boschker et al.)

- + Combined assessment of species and metabolic function
- Limited taxonomic resolution

2000: DNA-SIP (Radajewski et al.)

- + Combination of genetic methods allows even culture independent assessment of species
- Information on the potential but not the actual function
- Density gradient centrifugation provides only low resolution, results in high minimal labelling

2002: RNA-SIP (Manefield et al.)

- + Combination of genetic methods allows even culture independent assessment of species
- + Detection of transcription of functional genes
- Density gradient centrifugation enables only low resolution, results in high minimal labelling


2006: Nano-SIMS (Lechene et al.)

- + Single cell analysis
- + high dynamic range
- + Direct correlation between metabolomic activity and species identification in combination with FISH
- Due to tedious sample preparation only consortia of limited complexity can be analysed
- relatively costly instruments

2008: Protein-SIP (Jehmlich et al.)

- Labelled Protein can be extracted and analysed by gel-based or LC-MS-based approaches, mass spectrometry yields identification of protein (and species) and determination of metabolic activity
- Example of toluene degradation anoxic mixed culture
- + Mass spectrometry of peptides enables high sensitivity and dynamic range for incorporation
- + Direct correlation between function and metabolic activity
- Accuracy of annotation of proteins to a species depends on genomic information

# SIP-approaches in terms of taxonomic resolution, dynamic range and needed biomass





Further criteria can be : availability/cost of instrumentation  
time needed  
proof of functionality

# Workflow: Protein-based st<sub>a</sub>b<sub>i</sub>l<sub>e</sub> is<sub>t</sub>o<sub>p</sub>e pr<sub>o</sub>b<sub>i</sub>ng (Protein-SIP)


- „light“ substrate  
e.g.  $^{12}\text{C}_6$  benzene




- „heavy“ substrate  
e.g.  $^{13}\text{C}_6$  benzene




Protein extraction



tryptic digestion  
mass spectrometry



# Toluene degrading anoxic mixed culture for proof of principle




- anoxic degradation of  $^{13}\text{C}_7$  labeled toluene by *Azoarcus* EbN1  
In an artificial community with UFZ-1 culture (that does not grow on toluene but on gluconate)




Determination of  $^{13}\text{C}$  incorporation levels to link **functionality** and **identity** in the experimentally designed consortium

# EbN1 does not distinguish between $^{12}\text{C}$ and $^{13}\text{C}$ toluene



$^{13}\text{C}_7$  toluene: circles  
 $^{12}\text{C}_7$  toluene: triangles no  
substrate: dashed line


# UFZ-1 does not grow on toluene



gluconate: white diamonds  
toluene: black triangles

gluconate + toluene: black diamonds.

# Mixed culture of EbN1 and UFZ exhibits stable growth



$^{13}\text{C}_7$  toluene: white circles

$^{13}\text{C}_7$  toluene: black circles



gluconate: white diamonds

gluconate +  $^{12}\text{C}$  toluene: black diamonds





gluconate +  $^{13}\text{C}_7$ -toluene: grey triangles

No substrate: grey circles

# 2DE for pre-separation




# Molecular details in peptide mass fingerprints




CENTRE FOR  
ENVIRONMENTAL  
RESEARCH - UFZ

# Zoom in reveals altered isotopic envelopes





# Calculation of incorporation

## Measurement




## Computational simulation




➔ **98.7% incorporation**

# Molecular details on the MS/MS level





# $^{13}\text{C}$ incorporation takes place only into EbN1 and ...



detail view of  
one peptide

# ... not in UFZ-1 proteins



detail view of  
one peptide

# What about sensitivity and dynamic range of SIP approaches?

## → Sensitivity in respect to the amount of biomass

- so far most experiments are performed “off situ”
- DNA/RNA-SIP have the potential to be used “in situ” using backtraps with loaded substrates
- Protein-SIP might be possible with backtraps

## → Sensitivity in terms of detectable incorporation

- DNA/RNA-SIP are hampered by low resolution of density gradient centrifugation
- High resolution MS enables Protein-SIP to detect 2% incorporation of  $^{13}\text{C}$


## → Dynamic range

- DNA/RNA-SIP are mostly on/off
- Protein-SIP allows 2-100% range

# Composition of communities must be considered variable and dynamic

→ Who eats what, where and when?

→ Sensitivity and dynamic range must cope with variability of microbial communities



# Sensitivity and dynamic range of Protein-SIP

Aerobic growth of *P. putida* on  $^{13}\text{C}_6$  benzene with defined  $^{13}\text{C}$  content

$^{13}\text{C}$  10%




$^{13}\text{C}$  25%




$^{13}\text{C}$  50%



$^{13}\text{C}$  100%




# High resolution MS allows sensitivity down to 2% incorporation




With 10% incorporation the highest peak shifted 6 Da, indicating a sensitivity down to 2%

# Why do we need high sensitivity and dynamic range?

- Analysis of carbon flux in microbial communities
- Usage of incomplete labelled substrates reduces costs




1 mg Toluene  $^{13}\text{C}_6 = 435\text{€}$




1 mg Toluene  $^{13}\text{C}_1 = 12\text{-}45\text{€}$


- High sensitivity allows usage of specifically labelled positions



nonylphenol




MTBE



naphthalene

- Analysis of systems with high intrinsic dilution of substrate e.g. trophic networks



D-glucose



Uptake by  
bacteria




Uptake by  
protists


?

# Workflow: Protein-based stabile isotope probing (Protein-SIP)


- „light“ substrate  
e.g.  $^{12}\text{C}_6$  benzene






- „heavy“ substrate  
e.g.  $^{13}\text{C}_6$  benzene



Protein extraction



tryptic digestion  
mass spectrometry



# Half decimal place rule provides a comparison independent calculation of incorporation

→ half decimal place rule (HDPR):

mass values of tryptic peptides:

$$1,420.7 \rightarrow 14 : 2 = 7$$

$$1,630.8 \rightarrow 16 : 2 = 8$$

$$2,058.0 \rightarrow 20 : 2 = 10$$

$$2,389.2 \rightarrow 24 : 2 = 12$$


→  $^{12}\text{C}$  /  $^{13}\text{C}$  isotope mass differ by exactly 1.003355 amu

e.g.	1 C-atom	+0.003355 amu
	10 C-atoms	+0.03355 amu
	100 C-atoms	+0.3355 amu

→ accurate mass spectrometer is required


# Decimal slope shows regular distribution

- Helicobacter pylori tryptic peptides
- app. 80,000 data points
- each point: tryptic peptides fragments lengths varied between 2 and 40 amino acids



# Simulation points to a probable differentiation of $^{12}\text{C}$ and $^{13}\text{C}$ peptides


→ each  $^{12}\text{C}$  atom mass was substituted by their  $^{13}\text{C}$  atom mass to simulate a theoretical 100 atom % incorporation level




# Practical application using fully labeled $^{13}\text{C}$ *P. putida* peptides

→ calculated incorporation depends on the number of used peptide masses and the accuracy of MS measurement


$$^{13}\text{C}\% = \frac{\text{slope}(\text{exp}) - \text{slope}(^{12}\text{C})}{\text{slope}(^{13}\text{C}) - \text{slope}(^{12}\text{C})} \cdot 100\%$$



# The half decimal rule applies also to DNA/RNA



# 50 peptide masses are required for deviation < 5%



calculated incorporation depends on the number of used peptide masses and the accuracy of MS measurement

# HDPR yields comparison independent determination of incorporation

Aerobic growth of *P. putida* on  $^{13}\text{C}_6$  benzene with defined  $^{13}\text{C}$  content

$^{13}\text{C}$  10%











$^{13}\text{C}$  25%




$^{13}\text{C}$  50%




$^{13}\text{C}$  100%






# Workflow of Protein-SIP



# Protein-SIP can be extended to N and S



1-2 N/average amino acid  
5% of amino acids contain sulfur



# Conclusions

- Protein-SIP provides concomitantly information about functional proteins and metabolic activity of species
- Sensitivity in terms of biomass fits well with „*off situ*“ experiments, „*in situ*“ might be possible as well
- Protein-SIP, in contrast to DNA/RNA-SIP, features a dynamic range (2-100%) of incorporation of  $^{13}\text{C}$
- High resolution mass spectrometry and the Half decimal place rule allow a comparison independent and solely LC-MS based approach
- Protein-SIP can be extended to other isotopes ( $^{15}\text{N}$ ,  $^{36}\text{S}$ )
  
- Time resolved analysis of carbon flux in microbial communities becomes possible

# Example for importance of genomic information

*Cooperation with Prof. Suflita on Proteomics of alkane degraders*


*Desulfoglaeba alkanexedens*

Anaerobic sulfate reducer, degrades C<sub>6</sub>-C<sub>12</sub> *n*-alkanes completely to CO<sub>2</sub>


Isolated from a oily-wastewater storage facility

Cells were grown on butyrate and decane, both under sulfate-reducing conditions

# NCBI is not sufficient for environmentally important species



Growth with alkane:  
17 without and  
287 identifications  
with genome



Growth with butyrate:  
21 without and  
324 identifications  
with genome

# Further developments in Environmental Proteomics

- Improvement of isolation procedures
- Steady increase in genomic information of relevant species and communities
- Usage of proteins as Biomarkers for physiological properties in ecosystems
  
- Need for scientific exchange
- We plan to organize a meeting on Environmental Proteomics in 2011 in Leipzig to continue this event

# Acknowledgements

Department of Isotope Biogeochemistry  
Hans-Hermann Richnow  
Carsten Vogt  
Stephanie Hinke

Department of Environmental Microbiology  
Hauke Harms  
Susann Müller  
Thomas Hübschmann  
Roland Müller  
Ingo Fetzer

Department of Proteomics  
Nico Jehmlich  
Jana Seifert  
Martin Taubert  
Christine Schumann  
Michaela Öhler  
Yvonne Kullnick

University of Greifswald, Interfaculty Institute  
Genetics and Functional Genomics  
Frank Schmidt

MPI Institute for Infection Biology  
Jens Mattow

University of Oslo, Proteomics Group  
Bernd Thiede

University of Oklahoma  
Prof. Suflita  
Brendon Morris



This is a scanning electron micrograph (SEM) showing a highly detailed, porous, and interconnected network of biological or synthetic fibers. The structure is complex, with many small, rounded nodules and elongated, thin filaments. A semi-transparent grey rectangular box is superimposed over the center of the image, containing white text. At the bottom of the image, there is a black bar with white text providing technical details and a scale bar.

Hope to see you 2011 in Leipzig...

12C-Toluol + Nitrat Aussenseite 11.06.07 1kV 5mm

2 μm